

## CUSBSK-JP\_\_ Specification Sheet Basket Style



	P/N CUSBSK-JP		P/N CUSBSK-JPMAG	
Material	316 Stainless Steel		316 Stainless Steel with	
Waterial			Magnetic Stainless Steel	
Finish	Electropolished		Electropolished	
Retrieval System	Manual		Automated	
Overall Length	1.000 inches	25.4 mm	1.000 inches	25.4 mm
Inside Diameter	0.472 inches	12.0 mm	0.472 inches	12.0 mm
Wire Diameter	0.039 inches	1.0 mm	0.039 inches	1.0 mm
Pitch	0.125 inches	3.2 mm	0.125 inches	3.2 mm
Maximum Capsule Size	00		00	
Minimum Capsule Size	5		5	
Sinker Weight	2.15 grams		2.15 grams	

QLA Capsule Sinkers come in a wide range of sizes to fit all Dissolution applications, contact your QLA Sales Representative or visit our website at <a href="https://www.qla-llc.com">www.qla-llc.com</a> to view our current sinker product line.

2840 CLYMER AVE TELFORD, PA 18969 TEL: 908-685-7500 FAX: 908-526-4151 QLA-LLC.COM