


CAPWHT-XLS

Specification Sheet

Spiral Style



Material	316 Stainless Steel	
Finish	Electropolished	
Inside Length	1.150 inches	29.2 mm
Inside Diameter	0.465 inches	11.8 mm
Wire Diameter	0.043 inches	1.09 mm
Pitch	0.320 inches	8.1 mm
Maximum Width	0.690 inches	17.5 mm
Number of Coils	5.5	
Maximum Capsule Size	000	
Minimum Capsule Size	2	
Sinker Weight	2.15 grams	

QLA Capsule Sinkers come in a wide range of sizes to fit all Dissolution applications, contact your QLA Sales Representative or visit our website at www.qla-llc.com to view our current sinker product line.