

SNK010-VK Specification Sheet Basket Style



Material	316 Stainless Steel	
Finish	Electropolished	
Inside Length	0.620 inches	15.7 mm
Inside Diameter	0.795 inches	20.2 mm
Overall Height	0.800 inches	20.3 mm
Overall Diameter	1.000 inches	25.4 mm
Wire Diameter	0.025 inches	0.64 mm
Mesh Size	10 Mesh	
Maximum Capsule Size	1	
Minimum Capsule Size	5	
Sinker Weight	20.3 grams	
Agilent compatible P/N	12-3060	

QLA Capsule Sinkers come in a wide range of sizes to fit all Dissolution applications, contact your QLA Sales Representative or visit our website at www.qla-llc.com to view our current sinker product line.