


CAPWHT-NF

Specification Sheet

Spiral Style



Material	316 Stainless Steel	
Finish	PTFE Coated	
Inside Length	0.870 inches	22.0 mm
Inside Diameter	0.335 inches	8.5 mm
Wire Diameter (coated)	0.032 inches	0.81 mm
Pitch	0.158 inches	4.0 mm
Maximum Width	0.430 inches	11.0 mm
Number of Coils	4.2	
Maximum Capsule Size	1	
Minimum Capsule Size	5	
Sinker Weight	0.5 grams	

CAPWHT-NF is a spiral capsule sinker manufactured to meet the requirements of USP Nitrofurantoin Capsules Monograph.

QLA Capsule Sinkers come in a wide range of sizes to fit all Dissolution applications, contact your QLA Sales Representative or visit our website at www.qla-llc.com to view our current sinker product line.