

CAPWHT-S__ Series Specification Sheet O-Ring Style



	QLA P/N CAPWHT-S000		QLA P/N CAPWHT-S99	
Plate Body Material	316 Stainless Steel		316 Stainless Steel	
O-Ring Material	Viton		Viton	
Finish	Electropolished		Electropolished	
Inside Length	1.050 inches	26.6 mm	0.904 inches	23.0 mm
O-Ring Inside Diameter	0.301 inches	7.65 mm	0.301 inches	7.65 mm
Maximum Capsule Size	000		00	
Minimum Capsule Size	0		0	
Sinker Weight	1.15 grams		1.00 grams	

QLA Capsule Sinkers come in a wide range of sizes to fit all Dissolution applications, contact your QLA Sales Representative or visit our website at www.qla-llc.com to view our current sinker product line.

2840 CLYMER AVE TELFORD, PA 18969 TEL: 908-685-7500 FAX: 908-526-4151 QLA-LLC.COM