

## BSK008 Series Specification Sheet Basket Style



	BSK008-01		BSK008-01LG		BSK008-JP	
Material	316 Stainless Steel		316 Stainless Steel		316 Stainless Steel	
Finish	Electropolished		Electropolished		Electropolished	
Inside Length	1.06 Inches	26.9 mm	1.50 Inches	38.0 mm	0.90 Inches	22.9 mm
Inside Width	0.62 Inches	15.7 mm	0.62 Inches	15.7 mm	0.51 Inches	13.0 mm
Mesh Size	8 Mesh		8 Mesh		8 Mesh	
Maximum Capsule Size	000		000		0	
Minimum Capsule Size	5		5		5	
Lid Style	Clasp Type		Clasp Type		Clasp Type	
Sinker Weight	3.8 grams		6.1 grams		2.9 grams	

For dissolution testing of large capsules, tablets and suppositories, the BSK008 Series Sinkers are manufactured with precision tooling that produces consistent dimensions from sinker to sinker. This consistency of product will produce consistent Dissolution test results of your sample.

QLA Capsule Sinkers come in a wide range of sizes to fit all Dissolution applications, contact your QLA Sales Representative or visit our website at <a href="https://www.qla-llc.com">www.qla-llc.com</a> to view our current sinker product line.

2840 CLYMER AVE TELFORD, PA 18969 TEL: 908-685-7500 FAX: 908-526-4151 QLA-LLC.COM