


## CAPWHT-U\_\_ Specification Sheet 'U' Style



	QLA P/N CAPWHT-U01		QLA P/N CAPWHT-U2	
Material	316 Stainless Steel		316 Stainless Steel	
Finish	Electropolished		Electropolished	
Overall Height	1.024 inches	26.0 mm	0.984 inches	25.0 mm
Overall Width	0.750 inches	19.0 mm	0.692 inches	17.6 mm
Wire Diameter	0.031 inches	0.80 mm	0.031 inches	0.80 mm
Maximum Capsule Size	0		2	
Minimum Capsule Size	1		2	
Sinker Weight	0.34 grams		0.33 grams	

QLA Capsule Sinkers come in a wide range of sizes to fit all Dissolution applications, contact your QLA Sales Representative or visit our website at <a href="https://www.qla-llc.com">www.qla-llc.com</a> to view our current sinker product line.

2840 CLYMER AVE TELFORD, PA 18969 TEL: 908-685-7500 FAX: 908-526-4151 QLA-LLC.COM