


# CAPWST-31 Specification Sheet Spring Style



<b>Material</b>	316 Stainless Steel	
<b>Finish</b>	Electropolished	
<b>Inside Length</b>	1.220 inches	31.0 mm
<b>Inside Diameter</b>	0.430 inches	11.0 mm
<b>Wire Diameter</b>	0.031 inches	0.78 mm
<b>Number of Coils</b>	11	
<b>Maximum Capsule Size</b>	000	
<b>Minimum Capsule Size</b>	00	
<b>Sinker Weight</b>	3.25 grams	

**CAPWST-31 is compatible to Sotax P/N 7349**

QLA Capsule Sinkers come in a wide range of sizes to fit all Dissolution applications, contact your QLA Sales Representative or visit our website at [www.qla-llc.com](http://www.qla-llc.com) to view our current sinker product line.