


CAPWST-18 Specification Sheet Spring Style



Material	316 Stainless Steel	
Finish	Electropolished	
Inside Length	0.709 inches	18.0 mm
Inside Diameter	0.236 inches	6.0 mm
Wire Diameter	0.016 inches	0.42 mm
Number of Coils	8	
Maximum Capsule Size	3	
Minimum Capsule Size	4	
Sinker Weight	.75 grams	

CAPWST-18 is compatible to Sotax P/N 7346

QLA Capsule Sinkers come in a wide range of sizes to fit all Dissolution applications, contact your QLA Sales Representative or visit our website at www.qla-llc.com to view our current sinker product line.