


CAPWHT-4S

Specification Sheet

Spiral Style



Material	316 Stainless Steel	
Finish	Electropolished	
Inside Length	0.840 inches	21.3 mm
Inside Diameter	0.385 inches	9.8 mm
Wire Diameter	0.043 inches	1.09 mm
Pitch	0.280 inches	7.1 mm
Maximum Width	0.600 inches	15.2 mm
Number of Coils	8.5	
Maximum Capsule Size	00	
Minimum Capsule Size	3	
Sinker Weight	1.95 grams	

QLA Capsule Sinkers come in a wide range of sizes to fit all Dissolution applications, contact your QLA Sales Representative or visit our website at www.qla-llc.com to view our current sinker product line.