

CAPWHT-SV Specification Sheet Prong Style



Material	316 Stainless Steel	
Finish	Electropolished	
Outside Diameter (Large End)	0.365 inches	9.3 mm
Outside Length (Overall)	1.105 inches	28.0 mm
Inside Diameter (Small End)	0.210 inches	5.3 mm
Inside Diameter (Large End)	0.275 inches	7.0 mm
Capsule Sizes	1 through 3	
Sinker Weight	1.25 grams	

QLA Capsule Sinkers come in a wide range of sizes to fit all Dissolution applications, contact your QLA Sales Representative or visit our website at www.qla-llc.com to view our current sinker product line.