


CAPWHT-LG Specification Sheet Spiral Style



Material	Music Wire	
Finish	PTFE Coated	
Inside Length	1.100 inches	27.9 mm
Inside Diameter	0.410 inches	10.4 mm
Wire Diameter	0.043 inches	1.09 mm
Pitch	0.320 inches	8.1 mm
Maximum Width	0.620 inches	15.7 mm
Number of Coils	6.5	
Maximum Capsule Size	000	
Minimum Capsule Size	2	
Sinker Weight	2.0 grams	

QLA Capsule Sinkers come in a wide range of sizes to fit all Dissolution applications, contact your QLA Sales Representative or visit our website at www.qla-llc.com to view our current sinker product line.