


CAPWHT-0XS

Specification Sheet

Spiral Style



Material	Music Wire	
Finish	PTFE Coated	
Inside Length	0.953 inches	24.2 mm
Inside Diameter	0.575 inches	14.6 mm
Wire Diameter (coated)	0.047 inches	1.19 mm
Pitch	0.156 inches	4.0 mm
Maximum Width	0.790 inches	20.1 mm
Number of Coils	8.5	
Maximum Capsule Size	000	
Minimum Capsule Size	1	
Sinker Weight	3.29 grams	

QLA Capsule Sinkers come in a wide range of sizes to fit all Dissolution applications, contact your QLA Sales Representative or visit our website at www.qla-llc.com to view our current sinker product line.