

CUSBSK-ALF Specification Sheet Basket Style



Material	316L Stainless Steel	
Finish	Electropolished	
Overall Length	1.476 inches	37.5 mm
Outside Diameter	0.787 inches	20.0 mm
Cap Hole Diameter, Qty. 7	0.177 inches	4.5 mm
Bottom Hole Diameter, Qty. 7	0.177 inches	4.5 mm
Side Wall Hole Diameter, Qty. 60	0.197 inches	5.0 mm
Side Wall End Hole Diameter, Qty. 12	0.067 inches	1.7 mm
Sinker Weight	15.1 grams	

For dissolution testing of Alfuzosin Hydrochloride Extended Release Tablets, as specified in USP43-NF 38. The CUSBSK-ALF Tablet sinker is manufactured with precise dimensional hole placements to allow for exact replication of USP testing of the drug Alfuzosin Hydrochloride Extended Release tablets.

QLA Capsule Sinkers come in a wide range of sizes to fit all Dissolution applications, contact your QLA Sales Representative or visit our website at www.qla-llc.com to view our current sinker product line.