


# CAPWHT-02 Specification Sheet Spiral Style



<b>Material</b>	Music Wire	
<b>Finish</b>	PTFE Coated	
<b>Inside Length</b>	0.840 inches	21.3 mm
<b>Inside Diameter</b>	0.370 inches	9.4 mm
<b>Wire Diameter</b>	0.043 inches	1.09 mm
<b>Pitch</b>	0.156 inches	3.9 mm
<b>Maximum Width</b>	0.640 inches	16.2 mm
<b>Number of Coils</b>	8.5	
<b>Maximum Capsule Size</b>	00	
<b>Minimum Capsule Size</b>	4	
<b>Sinker Weight</b>	2.05 grams	

QLA Capsule Sinkers come in a wide range of sizes to fit all Dissolution applications, contact your QLA Sales Representative or visit our website at [www.qla-llc.com](http://www.qla-llc.com) to view our current sinker product line.