


# CAPWHT-LGMT Specification Sheet Spiral Style



<b>Material</b>	Music Wire	
<b>Finish</b>	Orange Non-Stick Epoxy Coated	
<b>Inside Length</b>	.944 inches	24.0 mm
<b>Inside Diameter</b>	0.485 inches	12.3 mm
<b>Wire Diameter</b>	0.043 inches	1.09 mm
<b>Pitch</b>	0.118 inches	3.00 mm
<b>Maximum Width</b>	0.650 inches	16.5 mm
<b>Number of Coils</b>	6.5	
<b>Maximum Capsule Size</b>	000	
<b>Minimum Capsule Size</b>	4	
<b>Sinker Weight</b>	1.9 grams	

QLA Capsule Sinkers come in a wide range of sizes to fit all Dissolution applications, contact your QLA Sales Representative or visit our website at [www.qla-llc.com](http://www.qla-llc.com) to view our current sinker product line.