


CAPWHT-1092__ Specification Sheet Spiral Style



	CAPWHT-10921&2		CAPWHT-10923&4	
Material	316 Stainless Steel		316 Stainless Steel	
Finish	Electropolished		Electropolished	
Inside Length	0.807 inches	20.5 mm	0.787 inches	20.0 mm
Inside Diameter	0.276 inches	7.0 mm	0.217 inches	5.5 mm
Wire Diameter	0.031 inches	0.80 mm	0.031 inches	0.80 mm
Pitch	0.197 inches	5.0 mm	0.193 inches	4.9 mm
Maximum Width	0.354 inches	9.0 mm	0.295 inches	7.5 mm
Number of Coils	4.1		4.0	
Maximum Capsule Size	1		3	
Minimum Capsule Size	2		4	
Sinker Weight	0.36 grams		0.32 grams	

CAPWHT-1092__ are spiral capsule sinkers manufactured to meet the requirements of USP, General Chapter, (1092) The Dissolution Procedure: Development and Validation, Table 2.

QLA Capsule Sinkers come in a wide range of sizes to fit all Dissolution applications, contact your QLA Sales Representative or visit our website at www.qla-llc.com to view our current sinker product line.