


CAPWSC-01S

Specification Sheet

Screw Cap Style



Material	316 Stainless Steel	
Finish	Electropolished	
Inside Length (Max)	0.765 inches	19.4 mm
Inside Diameter	0.281 inches	7.1 mm
Wire Diameter (Body)	0.060 inches	1.52 mm
Pitch	0.156 inches	4.0 mm
Maximum Width	0.411 inches	10.4 mm
Number of Coils	6	
Maximum Capsule Size	1	
Minimum Capsule Size	5	
Sinker Weight	3.25 grams	

QLA Capsule Sinkers come in a wide range of sizes to fit all Dissolution applications, contact your QLA Sales Representative or visit our website at www.qla-llc.com to view our current sinker product line.