


CAPWST-23 Specification Sheet Spring Style



Material	316 Stainless Steel	
Finish	Electropolished	
Inside Length	0.906 inches	23.0 mm
Inside Diameter	0.315 inches	8.0 mm
Wire Diameter	0.016 inches	0.42 mm
Number of Coils	9	
Maximum Capsule Size	0	
Minimum Capsule Size	1	
Sinker Weight	.95 grams	

CAPWST-23 is compatible to Sotax P/N 6881

QLA Capsule Sinkers come in a wide range of sizes to fit all Dissolution applications, contact your QLA Sales Representative or visit our website at www.qla-llc.com to view our current sinker product line.