



QLA SUPPOSITORY BASKETS



Certified Quality

QLA's Suppository Dissolution Baskets have been designed to meet the highest quality standards. These baskets are made from the high performance thermoplastic Polyoxymethylene (Delrin) and designed for improved structural integrity and maximum service life. Our Suppository Baskets are 100% inspected on an Optical Comparator and are guaranteed to provide 0.4 mm maximum total runout between the mounting bore and the lower rim. This feature is a critical component to ensure accurate and consistent results.

Ordering Information

<i>Instrument</i>	<i>QLA Part Number</i>
Agilent/VanKel/Varian	BSKSUP-01
Caleva	BSKSUP-CA
Distek Standard Push-On	BSKSUP-DK
Distek Evolution	BSKSUP-DKC
Erweka	BSKSUP-EW
PharmaTest	BSKSUP-PT
Sotax	BSKSUP-ST
Teledyne Hanson Research	BSKSUP-HR

Contact your QLA Sales Representative for any special requirements.

YOUR DISSOLUTION SPECIALISTS