

# SNK008-VK Specification Sheet Basket Style



<b>Material</b>	316 Stainless Steel	
<b>Finish</b>	Electropolished	
<b>Inside Length</b>	0.620 inches	15.7 mm
<b>Inside Diameter</b>	0.795 inches	20.2 mm
<b>Overall Height</b>	0.800 inches	20.3 mm
<b>Overall Diameter</b>	1.000 inches	25.4 mm
<b>Wire Diameter</b>	0.028 inches	0.71 mm
<b>Mesh Size</b>	8 x 8	
<b>Maximum Capsule Size</b>	1	
<b>Minimum Capsule Size</b>	5	
<b>Sinker Weight</b>	20.3 grams	
<b>Agilent compatible P/N</b>	12-3062	

QLA Capsule Sinkers come in a wide range of sizes to fit all Dissolution applications, contact your QLA Sales Representative or visit our website at [www.qla-llc.com](http://www.qla-llc.com) to view our current sinker product line.