


CAPWHT-XXLS

Specification Sheet

Spiral Style



Material	316 Stainless Steel	
Finish	Electropolished	
Inside Length	1.600 inches	40.6 mm
Inside Diameter	0.512 inches	13.0 mm
Wire Diameter	0.043 inches	1.09 mm
Pitch	0.310 inches	7.9 mm
Maximum Width	0.730 inches	18.5 mm
Number of Coils	6	
Maximum Capsule Size	Special	
Minimum Capsule Size	0	
Sinker Weight	3.0 grams	

QLA Capsule Sinkers come in a wide range of sizes to fit all Dissolution applications, contact your QLA Sales Representative or visit our website at www.qla-llc.com to view our current sinker product line.