


Distek 4300 Syringe Pump Tubing Replacement Kits




TUBSP-DK4300



TUBSP-DK4300-6V



TUBSP-DK4300-6R



TUBSP-DK4300-6C

All kits come complete with integrated silicone pinch-valve tubing sections and color-coded fittings for easy installation.

TUBSP-DK4300 - Tubing Set for 6 Position Distek Evolution 4300 Syringe Pump. Includes 18 lines with fittings for vessel, auto-sampler and replacement sections.

TUBSP-DK4300-6V- Tubing Set for 6 Position Distek Evolution 4300 Syringe Pump. Includes 6 lines with fittings for vessel section.

TUBSP-DK4300-6R- Tubing Set for 6 Position Distek Evolution 4300 Syringe Pump. Includes 6 lines with fittings for replacement section.

TUBSP-DK4300-6C - Tubing Set for 6 Position Distek Evolution 4300 Syringe Pump. Includes 6 lines with fittings for collector/auto-sampler section.

Wetted Materials

Component	Material
Fitting (all colors)	Delrin™
Ferrule (Blue)	PTFE
Sample Tubing	PTFE
Pinch Valve tubing	Silicone

For more information please contact your QLA Representative or Distributor.