

# New Product Releases



## **3<sup>RD</sup> QTR 2023**



QUALITY LAB ACCESSORIES

## 970-02-050

DDAPT-CH Replacement Arm for  
use with APEX Vessels





## BSK040-ST23FIN

40 mesh Clip Style Basket with 3 fins for Sotax,  
Rounded Rim Sotax style, 316 SS, Serialized  
(15 degree axial fins)

## BSK100-ST2

100 Mesh Clip Style Basket for Sotax,  
Rounded Rim Sotax style, 316 SS,  
Serialized





## BSKPTEW-40

40 Mesh (USP) Clip Style Basket for  
Erweka, PTFE Coated over 316 SS,  
Serialized

## CUSBSK-20/200EL

200 Mesh Push-On Style Basket with  
20 Mesh Stability Lining for Electrolab,  
316 SS, Serialized





## DISTUB-DK06

Glass Replacement Tubes for 6 Tube  
Disintegration Assembly for Distek Testers,  
25mm OD x 77.5mm L (Set/6)

## FIL001-EL

1 Micron Porous Filters, UHMW Polyethylene, 3mm ID, Electrolab compatible (Jar/1000)





## FIL001-EL-a

1 Micron Porous Filters, UHMW Polyethylene, 3mm ID, Electrolab compatible (Pack/100)

## FIL001-LG

1 Micron Porous Filters, UHMW  
Polyethylene, Logan compatible  
(Jar/1000)





## FIL001-LG-a

1 Micron Porous Filters,  
UHMW Polyethylene, Logan  
compatible (Pack/100)

## FIL005-EL

5 Micron Porous Filters, UHMW  
Polyethylene, 3mm ID, Electrolab  
compatible (Jar/1000)





## FIL005-EL-a

5 Micron Porous Filters, UHMW  
Polyethylene, 3mm ID, Electrolab  
compatible (Pack/100)

## FIL005-LG

5 Micron Porous Filters, UHMW  
Polyethylene, Logan compatible  
(Jar/1000)





## FIL005-LG-a

5 Micron Porous Filters, UHMW Polyethylene,  
Logan compatible (Pack/100)



## FILTOOL-HR

Filter ID Adjustment Tool for Hanson  
1/8" ID Filters (if needed when used  
with PEEK 1/8" OD cannula)




## KNOB-HR

Vessel Clamp Knob for Hanson  
SR8-Plus baths (set of 2)

## PRNS10-ST62

900mL 316 SS Cannula without  
Filter Housing (.062"/1.6mm  
tubing, 10"/254mm length) for  
Sotax AT 7/ AT7smart/ AT Xtend





**Please speak with your Sales Rep regarding  
availability of the above.**



QUALITY LAB ACCESSORIES