

New Product Releases


2nd QTR 2024



QUALITY LAB ACCESSORIES

CANS09-SH

9" (228.6mm) Straight 316
SS 1/8" (3.2mm) OD Cannula
with SS Luer Lock





CAPWHT-10920

Spiral Capsule Sinker for Capsule Size 0 per
USP <1092>, 316 SS, Approx .93" (23.7mm) L x
.33" (8.5mm) W capacity, 4 coils (.032"/20
gauge wire diameter (Initial wire length=12cm))

DEPSET-PK2L

Depth Set Tool for Apex 2L vessels,
Serialized with 1 year calibration





FIL020-VK

20 Micron Porous Filters, UHMW Polyethylene,
Direct Fit-No Tygon Tubing Needed, Agilent/
VanKel/Varian compatible (Jar/1000)

FIL020-VK-a

20 Micron Porous Filters, UHMW Polyethylene,
Direct Fit-No Tygon Tubing Needed, Agilent/
VanKel/Varian compatible (Pack/100)





GLA900-STXGL

1000mL Clear Glass Vessel with Graduation
Lines for Sotax Xtend, Serialized

NPD013-EW

13.5" (343mm) Paddle with 10mm OD for Erweka, PTFE coated over 316 SS, Serialized





TRANDEHLD-2.5

Transdermal Patch Holder, 2.5 cm sq, 316 SS

TRANDEHLD-5

Transdermal Patch Holder, 5 cm sq, 316 SS




TUBSP-DK4300-8RBN

Tubing Set for 8 Position Distek Evolution 4300 Syringe Pump. Includes 8 lines with fittings for Recycle Bin and Needle Assembly, 53" lengths (Set/8)

TUBSP- DK4300-8RCD

Tubing Set for 8 Position Distek Evolution 4300 Syringe Pump. Includes 8 lines with fittings from Recycle to Diverter, 24" lengths (Set/8)





**Please speak with your Sales Rep regarding
availability of the above.**



QUALITY LAB ACCESSORIES